

## COUPLAND SUCTION TUBE INSERT

Ø 2.3mm; 95mm



**SKU:** 085-010-23

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>AUTHORS</b>	<a href="#">COUPLAND</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>AUTHORS</b>	<a href="#">COUPLAND</a>