

**COUPLAND SUCTION TUBE INSERT**

Ø 1.5mm; 95mm

**SKU:** 085-010-15**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)**ZUSÄTZLICHE INFORMATIONEN**

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>KLS Martin</b>	18-575-05
<b>AUTHORS</b>	<a href="#">COUPLAND</a>

**ZUSÄTZLICHE INFORMATIONEN**

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>KLS Martin</b>	18-575-05
<b>AUTHORS</b>	<a href="#">COUPLAND</a>