

## CRILE NEEDLE HOLDER

straight; serrated; without catch

15.0cm (6")

**SKU:** 050-338-15

**Categories:** [Open Surgery](#), [Suture](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM503
<b>KLS Martin</b>	20-500-15
<b>AUTHORS</b>	<a href="#">CRILE</a>
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM503
<b>KLS Martin</b>	20-500-15
<b>AUTHORS</b>	<a href="#">CRILE</a>

## SPECIFICATIONS

[serrated](#), [straight](#)