

CW-CASPAR DISC RONGEUR

angled up; 4 x 14mm; toothed

18cm (7")

SKU: 240-515-04-CW

Categories: [Neurosurgery](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Laminectomy Rongeurs](#)
AUTHORS [CASPAR](#), [CLEANWAVE](#)
SPECIFICATIONS [angled up](#), [toothed](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Laminectomy Rongeurs](#)
AUTHORS [CASPAR](#), [CLEANWAVE](#)
SPECIFICATIONS [angled up](#), [toothed](#)