

## CW-CASPAR DISC RONGEUR

angled down; 2 x 12mm; toothed

18.5cm (7.4")

**SKU:** 240-518-02-CW

**Categories:** [Neurosurgery](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Laminectomy Rongeurs</a>
<b>AUTHORS</b>	<a href="#">CASPAR</a> , <a href="#">CLEANWAVE</a>
<b>SPECIFICATIONS</b>	<a href="#">angled down</a> , <a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Laminectomy Rongeurs</a>
<b>AUTHORS</b>	<a href="#">CASPAR</a> , <a href="#">CLEANWAVE</a>
<b>SPECIFICATIONS</b>	<a href="#">angled down</a> , <a href="#">toothed</a>