

CW-CASPAR DISC RONGEUR

angled down; 3 x 12mm; toothed

18.5cm (7.4")

SKU: 240-518-03-CW

Categories: [Neurosurgery](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Laminectomy Rongeurs](#)
AUTHORS [CASPAR](#), [CLEANWAVE](#)
SPECIFICATIONS [angled down](#), [toothed](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Laminectomy Rongeurs](#)
AUTHORS [CASPAR](#), [CLEANWAVE](#)
SPECIFICATIONS [angled down](#), [toothed](#)