

CW-CASPAR DISC RONGEUR

straight; 2 x 12mm; toothed

18cm (7")

SKU: 240-511-02-CW

Categories: [Open Surgery](#), [Neurosurgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Rongeurs
AUTHORS	CASPAR , CLEANWAVE
SPECIFICATIONS	straight , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Rongeurs
AUTHORS	CASPAR , CLEANWAVE
SPECIFICATIONS	straight , toothed