

DALE MULTI PURPOSE VASCULAR CLAMP

angled; atraumatic

18.0cm (7.0")




SKU: 142-025-18

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
KLS Martin	24-340-18
AUTHORS	DALE
Medicon	55.52.51
SPECIFICATIONS	angled



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
KLS Martin	24-340-18
AUTHORS	DALE
Medicon	55.52.51

SPECIFICATIONS

[angled](#)