

## DANDY-BABY HAEMOSTATIC FORCEPES

curved; serrated

12.0cm (4.8")



**SKU:** 040-341-12

**Categories:** [Open Surgery](#), [Hemostats](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>KLS Martin</b>	13-351-12
<b>AUTHORS</b>	<a href="#">BABY</a> , <a href="#">DANDY</a>
<b>Medicon</b>	16.11.83
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>KLS Martin</b>	13-351-12
<b>AUTHORS</b>	<a href="#">BABY</a> , <a href="#">DANDY</a>
<b>Medicon</b>	16.11.83

## SPECIFICATIONS

[curved](#), [serrated](#)