

DARDIK MULTI PURPOSE VASCULAR CLAMP

90° angled; atraumatic

13.5cm (5.3")




SKU: 142-022-90

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---------------------------------|
| UNIT | pcs. |
| TYPE | Vascular Clamps |
| Aesculap | FC314 |
| KLS Martin | 24-338-90 |
| AUTHORS | DARDIK |
| SPECIFICATIONS | angled |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|---------------------------------|
| UNIT | pcs. |
| TYPE | Vascular Clamps |
| Aesculap | FC314 |
| KLS Martin | 24-338-90 |
| AUTHORS | DARDIK |

SPECIFICATIONS

[angled](#)