

DARDIK MULTI PURPOSE VASCULAR CLAMP

15° angled; atraumatic

15.5cm (6.1")



SKU: 142-020-15

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
Aesculap	FC301
KLS Martin	24-336-15
AUTHORS	DARDIK
SPECIFICATIONS	angled



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
Aesculap	FC301
KLS Martin	24-336-15
AUTHORS	DARDIK

SPECIFICATIONS

[angled](#)