

DARDIK MULTI PURPOSE VASCULAR CLAMP

30° angled; atraumatic

15.0cm (5.9")



SKU: 142-020-30

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
Aesculap	FC302
KLS Martin	24-336-30
AUTHORS	DARDIK
SPECIFICATIONS	angled



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
Aesculap	FC302
KLS Martin	24-336-30
AUTHORS	DARDIK

SPECIFICATIONS

[angled](#)