

DARDIK MULTI PURPOSE VASCULAR CLAMP


90° angled; atraumatic

13.5cm (5.3")



SKU: 142-020-90

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
Aesculap	FC304
KLS Martin	24-336-90
AUTHORS	DARDIK
SPECIFICATIONS	angled

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
Aesculap	FC304
KLS Martin	24-336-90
AUTHORS	DARDIK

SPECIFICATIONS

[angled](#)