

DARDIK MULTI PURPOSE VASCULAR CLAMP

straight; atraumatic

17.0cm (6.6")



SKU: 142-022-00

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
Aesculap	FC310
KLS Martin	24-338-00
AUTHORS	DARDIK
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
Aesculap	FC310
KLS Martin	24-338-00
AUTHORS	DARDIK

SPECIFICATIONS

[straight](#)