

DARDIK MULTI PURPOSE VASCULAR CLAMP

30° angled; atraumatic

16.5cm (6.4")



SKU: 142-022-30

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
Aesculap	FC312
KLS Martin	24-338-30
AUTHORS	DARDIK
SPECIFICATIONS	angled



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
Aesculap	FC312
KLS Martin	24-338-30
AUTHORS	DARDIK

SPECIFICATIONS

[angled](#)