

DARDIK MULTI PURPOSE VASCULAR CLAMP

45° angled; atraumatic

17.0cm (6.6")




SKU: 142-022-45

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
Aesculap	FC313
KLS Martin	24-338-45
AUTHORS	DARDIK
SPECIFICATIONS	angled



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
Aesculap	FC313
KLS Martin	24-338-45
AUTHORS	DARDIK

SPECIFICATIONS

[angled](#)