

DARDIK MULTI PURPOSE VASCULAR CLAMP

straight; atraumatic

15.0cm (5.9")



SKU: 142-020-00

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
Aesculap	FC300
KLS Martin	24-336-00
AUTHORS	DARDIK
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
Aesculap	FC300
KLS Martin	24-336-00
AUTHORS	DARDIK

SPECIFICATIONS

[straight](#)