

DE BAKEY DIETHRICH BULLDOG CLAMP

straight; length: 50mm; jaw 14mm

closing force ca. 180g

SKU: 026-022-02

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	FB364
KLS Martin	13-122-14
AUTHORS	DE BAKEY , DIETRICH
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	FB364
KLS Martin	13-122-14
AUTHORS	DE BAKEY , DIETRICH

SPECIFICATIONS

[straight](#)