

DE BAKEY DIETHRICH BULLDOG CLAMP

straight; length: 60mm; jaw 20mm

closing force ca. 180g

SKU: 026-022-03

Categories: [Open Surgery](#), [Hemostats](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	FB366
KLS Martin	13-122-20
AUTHORS	DE BAKEY , DIETRICH
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	FB366
KLS Martin	13-122-20
AUTHORS	DE BAKEY , DIETRICH

SPECIFICATIONS

[straight](#)