

## DE BAKEY DIETHRICH BULLDOG CLAMP

angled; length: 45mm; jaw 14mm

closing force ca. 180g


**SKU:** 026-023-02

**Categories:** [Open Surgery](#), [Hemostats](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bulldog Clamps</a>
<b>Aesculap</b>	FB365
<b>KLS Martin</b>	13-123-14
<b>AUTHORS</b>	<a href="#">DE BAKEY</a> , <a href="#">DIETRICH</a>
<b>SPECIFICATIONS</b>	<a href="#">angled</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bulldog Clamps</a>
<b>Aesculap</b>	FB365
<b>KLS Martin</b>	13-123-14
<b>AUTHORS</b>	<a href="#">DE BAKEY</a> , <a href="#">DIETRICH</a>

## SPECIFICATIONS

[angled](#)