

DE BAKEY-SATINSKY BULLDOG CLAMP

angled; atraumatic

10.0cm (4")


SKU: 026-191-10

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	FB358
KLS Martin	24-291-10
AUTHORS	DE BAKEY , SATINSKY
SPECIFICATIONS	angled , atraumatic



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	FB358
KLS Martin	24-291-10
AUTHORS	DE BAKEY , SATINSKY

SPECIFICATIONS

[angled](#), [atraumatic](#)