

## DE BAKEY SUCTION TUBE

removable suction tip; Luer connector

inner Ø 3.0mm; outer Ø 6.0mm; 27.0cm

**SKU:** 085-008-27

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>Aesculap</b>	GF946
<b>AUTHORS</b>	<a href="#">DE BAKEY</a>
<b>Medicon</b>	09.85.07

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>Aesculap</b>	GF946
<b>AUTHORS</b>	<a href="#">DE BAKEY</a>
<b>Medicon</b>	09.85.07