

DERF-CONVERSE NEEDLE HOLDER

straight; longitudinal serrated

12.5cm (5")




SKU: 050-020-12

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM200
KLS Martin	20-328-12
AUTHORS	CONVERSE , DERF
Medicon	10.12.82
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM200
KLS Martin	20-328-12

AUTHORS

[CONVERSE](#), [DERF](#)

Medicon

10.12.82

SPECIFICATIONS

[serrated](#), [straight](#)