

## DIAMOND NITINOL LITHOTRIPTERBASKET

20 x 40mm; 8 wires

2.6mm; 230cm



**SKU:** PK-1LDXX2620230AK

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Stonebaskets for Retrieval & Lithotripsy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Stonebaskets for Retrieval & Lithotripsy](#)