

DIEFFENBACH BULLDOG CLAMP

straight; length: 38mm

SKU: 026-014-03

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)


Aesculap BH020

KLS Martin 13-114-03

AUTHORS [DIEFFENBACH](#)

Medicon 15.10.38

SPECIFICATIONS [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

Aesculap BH020

KLS Martin 13-114-03

AUTHORS [DIEFFENBACH](#)

Medicon 15.10.38

SPECIFICATIONS

[straight](#)