

## DIEFFENBACH BULLDOG CLAMP

straight; length: 38mm

**SKU:** 026-014-03

**Categories:** [Open Surgery](#), [Hemostats](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Bulldog Clamps](#)


**Aesculap** BH020

**KLS Martin** 13-114-03

**AUTHORS** [DIEFFENBACH](#)

**Medicon** 15.10.38

**SPECIFICATIONS** [straight](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Bulldog Clamps](#)

**Aesculap** BH020

**KLS Martin** 13-114-03

**AUTHORS** [DIEFFENBACH](#)

**Medicon** 15.10.38

## SPECIFICATIONS

[straight](#)