

DIETHRICH BULLDOG CLAMP

straight; length: 60mm; jaw 20mm

closing force ca. 80g

SKU: 026-020-03

Categories: [Open Surgery](#), [Hemostats](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
KLS Martin	13-120-12
AUTHORS	DIETHRICH
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
KLS Martin	13-120-12
AUTHORS	DIETHRICH
SPECIFICATIONS	straight