

## DIETHRICH BULLDOG CLAMP

straight; length: 60mm; jaw 20mm

closing force ca. 80g

**SKU:** 026-020-03

**Categories:** [Open Surgery](#), [Hemostats](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bulldog Clamps</a>
<b>KLS Martin</b>	13-120-12
<b>AUTHORS</b>	<a href="#">DIETRICH</a>
<b>SPECIFICATIONS</b>	<a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bulldog Clamps</a>
<b>KLS Martin</b>	13-120-12
<b>AUTHORS</b>	<a href="#">DIETRICH</a>
<b>SPECIFICATIONS</b>	<a href="#">straight</a>