

DIETHRICH BULLDOG CLAMP

straight; length: 60mm; jaw 20mm

closing force ca. 80g

SKU: 026-020-03

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

KLS Martin 13-120-12

AUTHORS [DIETHRICH](#)

SPECIFICATIONS [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

KLS Martin 13-120-12

AUTHORS [DIETHRICH](#)

SPECIFICATIONS [straight](#)