

DIETHRICH BULLDOG CLAMP

angled; length: 45mm; jaw 12mm

closing force ca. 80g

SKU: 026-021-02

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

Aesculap BH034

KLS Martin 13-121-10

AUTHORS [DIETRICH](#)

SPECIFICATIONS [angled](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

Aesculap BH034

KLS Martin 13-121-10

AUTHORS [DIETRICH](#)

SPECIFICATIONS

[angled](#)