

DIETHRICH BULLDOG CLAMP

angled; length: 50mm; jaw 20mm

closing force ca. 80g


SKU: 026-021-03

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
KLS Martin	13-121-12
AUTHORS	DIETRICH
SPECIFICATIONS	angled



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
KLS Martin	13-121-12
AUTHORS	DIETRICH
SPECIFICATIONS	angled