

DIETHRICH BULLDOG CLAMP

straight; length 45mm; jaw 8mm

closing force ca. 80g

SKU: 026-020-01

Categories: [Open Surgery](#), [Hemostats](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
KLS Martin	13-120-08
AUTHORS	DIETHRICH
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
KLS Martin	13-120-08
AUTHORS	DIETHRICH
SPECIFICATIONS	straight