

DOYEN INTESTINAL CLAMP

straight; diagonally serrated

23.0cm (9.2")



SKU: 160-210-23

Categories: [Open Surgery](#), [Stomach](#), [Intestines](#), [Rectum](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Stomach and Intestinal Clamp Forceps
Aesculap	EA126
KLS Martin	30-274-23
AUTHORS	DOYEN
Medicon	50.16.23
SPECIFICATIONS	diagonally serrated , straight

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Stomach and Intestinal Clamp Forceps
Aesculap	EA126
KLS Martin	30-274-23

AUTHORS

[DOYEN](#)

Medicon

50.16.23

SPECIFICATIONS

[diagonally serrated](#), [straight](#)