

## DOYEN INTESTINAL CLAMP

straight; diagonally serrated

23.0cm (9.2")



**SKU:** 160-210-23

**Categories:** [Open Surgery](#), [Stomach](#), [Intestines](#), [Rectum](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Stomach and Intestinal Clamp Forceps</a>
<b>Aesculap</b>	EA126
<b>KLS Martin</b>	30-274-23
<b>AUTHORS</b>	<a href="#">DOYEN</a>
<b>Medicon</b>	50.16.23
<b>SPECIFICATIONS</b>	<a href="#">diagonally serrated</a> , <a href="#">straight</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Stomach and Intestinal Clamp Forceps</a>
<b>Aesculap</b>	EA126
<b>KLS Martin</b>	30-274-23

## AUTHORS

[DOYEN](#)

## Medicon

50.16.23

## SPECIFICATIONS

[diagonally serrated](#), [straight](#)