

## DOYEN MALLET

solid; head 35mm; 240g

21.0cm (8.4")


**SKU:** 191-212-35

**Categories:** [Open Surgery](#), [Orthopaedics](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Mallets</a>
<b>KLS Martin</b>	23-204-03
<b>AUTHORS</b>	<a href="#">DOYEN</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Mallets</a>
<b>KLS Martin</b>	23-204-03
<b>AUTHORS</b>	<a href="#">DOYEN</a>