

## DOYEN SPHERICAL BURR

short Borchardt shaft

Ø 16.0mm



**SKU:** 190-291-16

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Hand Drills</a>
<b>KLS Martin</b>	22-212-16
<b>AUTHORS</b>	<a href="#">DOYEN</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Hand Drills</a>
<b>KLS Martin</b>	22-212-16
<b>AUTHORS</b>	<a href="#">DOYEN</a>