

DOYEN SPHERICAL BURR

short Borchardt shaft

Ø 4.0mm



SKU: 190-291-04

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hand Drills
KLS Martin	22-212-04
AUTHORS	DOYEN

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hand Drills
KLS Martin	22-212-04
AUTHORS	DOYEN