

## DOYEN SPHERICAL BURR

long Borchardt shaft

Ø 4.0mm



**SKU:** 190-290-04

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Hand Drills</a>
<b>KLS Martin</b>	22-210-04
<b>AUTHORS</b>	<a href="#">DOYEN</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Hand Drills</a>
<b>KLS Martin</b>	22-210-04
<b>AUTHORS</b>	<a href="#">DOYEN</a>