

## DOYEN TONGUE SPATULA

double ended; 23x30mm

13.0cm (5.2")

**SKU:** 114-463-13

**Categories:** [Open Surgery](#), [Tonsils](#), [Tracheotomy](#), [Tongue](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Cheek and Lip Retractors, Tongue Depressors</a>
<b>KLS Martin</b>	38-375-13
<b>AUTHORS</b>	<a href="#">DOYEN</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Cheek and Lip Retractors, Tongue Depressors</a>
<b>KLS Martin</b>	38-375-13
<b>AUTHORS</b>	<a href="#">DOYEN</a>