

DUPLEX-BABY STETHOSCOPES

Ø 36mm



SKU: 100-322-00

Categories: [Open Surgery](#), [Diagnostics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Stethoscopes
Aesculap	AC150
KLS Martin	17-260-00
AUTHORS	DUPLEX

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Stethoscopes
Aesculap	AC150
KLS Martin	17-260-00
AUTHORS	DUPLEX