

## DUPLEX-BABY STETHOSCOPES

Ø 36mm



**SKU:** 100-322-00

**Categories:** [Open Surgery](#), [Diagnostics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Stethoscopes</a>
<b>Aesculap</b>	AC150
<b>KLS Martin</b>	17-260-00
<b>AUTHORS</b>	<a href="#">DUPLEX</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Stethoscopes</a>
<b>Aesculap</b>	AC150
<b>KLS Martin</b>	17-260-00
<b>AUTHORS</b>	<a href="#">DUPLEX</a>