

## EAR NEEDLE

angled; 90° upwards; 1.0mm

16.0cm (6.4")

**SKU:** 111-648-10

**Categories:** [Open Surgery](#), [Otology](#)



## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needles</a>
<b>KLS Martin</b>	36-025-03



## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needles</a>
<b>KLS Martin</b>	36-025-03