

ENGLISH MODEL TOWEL FORCEPS

sharp
 9.0cm (3.6")



SKU: 030-109-09
Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Towel Forceps
Aesculap	BF408
KLS Martin	13-877-09
AUTHORS	ENGLISCHES MODELL
Medicon	17.50.28
SPECIFICATIONS	curved , sharp

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Towel Forceps
Aesculap	BF408
KLS Martin	13-877-09

AUTHORS

[ENGLISCHES MODELL](#)

Medicon

17.50.28

SPECIFICATIONS

[curved](#), [sharp](#)