

## ERCPC-CATHETER; BOTTLE TOP SHAPED TIP

tube 1.8mm; 210cm (single use)

able for guide wire 0.035"



**SKU:** PK-3CXEW1812210M-S

**Category:** [Endoscopy](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)