

ERCP-CATHETER; BOTTLE TOP SHAPED TIP

tube: 1.8mm/210cm

able for guide wire 0.035"

SKU: PK-3CXPk1810210M

Category: [Endoscopy](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [ERCP-Catheter, Papillotomes](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [ERCP-Catheter, Papillotomes](#)