

ERCP-CATHETER; BOTTLE TOP SHAPED TIP

able for guide wire 0.035"

.8mm, length 250cm single use



SKU: PK-3CXEW1810250M-S

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Instruments for Double-Balloon-Enteroscopia](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Instruments for Double-Balloon-Enteroscopia](#)