

## FERGUSSON (POPPEN) SUCTION TUBE

11 Fr.; 45° angled

14.5cm (5.8") working length



**SKU:** 085-034-11

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>AUTHORS</b>	<a href="#">FERGUSSON, POPPEN</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>AUTHORS</b>	<a href="#">FERGUSSON, POPPEN</a>