

## FERROZELL LIGATURE CONDUCTOR

straight

24.0cm (9.6")

**SKU:** 056-520-24

**Categories:** [Open Surgery](#), [Suture](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Ligature Conductors</a>
<b>Aesculap</b>	BM938
<b>KLS Martin</b>	20-870-24
<b>AUTHORS</b>	<a href="#">FERROZELL</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Ligature Conductors</a>
<b>Aesculap</b>	BM938
<b>KLS Martin</b>	20-870-24
<b>AUTHORS</b>	<a href="#">FERROZELL</a>