

## FINOCHIETTO NEEDLE HOLDER

angled; serrated; grooved jaw

27.0cm (10.8")


**SKU:** 050-211-27

**Categories:** [Open Surgery](#), [Suture](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM277
<b>KLS Martin</b>	20-415-27
<b>AUTHORS</b>	<a href="#">FINOCHIETTO</a>
<b>Medicon</b>	10.43.27
<b>SPECIFICATIONS</b>	<a href="#">angled</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM277
<b>KLS Martin</b>	20-415-27

## AUTHORS

[FINOCHIETTO](#)

## Medicon

10.43.27

## SPECIFICATIONS

[angled](#), [serrated](#)