

## FINSEN SELFRETAINING RETRACTOR

blunt; 2X3 teeth

5.0cm (2")


**SKU:** 070-611-05

**Categories:** [Open Surgery](#), [Retractors](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Selfretaining Retractors</a>
<b>KLS Martin</b>	15-701-05
<b>AUTHORS</b>	<a href="#">FINSEN</a>
<b>Medicon</b>	24.04.06
<b>SPECIFICATIONS</b>	<a href="#">blunt</a> , <a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Selfretaining Retractors</a>
<b>KLS Martin</b>	15-701-05
<b>AUTHORS</b>	<a href="#">FINSEN</a>
<b>Medicon</b>	24.04.06

## SPECIFICATIONS

[blunt](#), [toothed](#)